

LPS – Lab Productivity Systems

Xsample 340

Automation... what for?

QUALITY AND PRODUCTIVITY MATTERS

Why do automation in labs?

► Reduce costs per measurement...


- Lab staff can do more important things than sample handling
- Automated measurement even of high-viscous samples



► Reduce handling errors and increase uptime of instruments...

- User independent filling of master instrument
- Automated cleaning with pre-defined workflows





A first look at the instrument

ONE DEVICE – DIFFERENT SYRINGES

- ▶ Automated filling with syringes
- ▶ Active cleaning with 2 rinsing agents
- ▶ Force sensor to adjust filling speed
- ▶ User-independent repeatability
- ▶ Time savings at high-viscous samples due to automation

WELL DONE BUT ...



► Adapters for your favorite syringes

- 2 mL Braun syringes
- 5 mL Braun syringes
- 10 mL Braun syringes
- 5 mL ILS glass syringes

... by a simple exchange of an adapter

► Outstanding robustness

... trouble free operation under harsh conditions



Where can these features be used.

APPLICATIONS

- ▶ Syringe is the preferred consumable
- ▶ Automatic rinsing and drying is needed (cleaning module)
- ▶ Less samples per day need to be measured
- ▶ High-viscous samples are filled into the measuring system

What can the Xsample 340 be used for.

PETROLEUM INDUSTRY

FOOD INDUSTRY

PHARMACEUTICAL INDUSTRY

Special challenges for the sampler unit...

- ▶ High-viscous samples
- ▶ Easy to clean
- ▶ Volatile samples
- ▶ Harsh conditions
- ▶ Easy retrofit of systems

Our solution...



Do I need to set the filling speed to a certain value to get high-viscous sample filled?

Force sensor controlled filling:

- ▶ No, because the system automatically adjusts the speed to the maximum allowed level. The force-controlled movement of the pusher prevents excess pressure in the system and ensures smooth filling.




Easy to clean...

Do I need to get in contact with cleaning agents after each measurement?

Easy automated cleaning:

- ▶ The whole system is cleaned with up to two different cleaning agents, then dried automatically. Once the cleaning agents are filled and in place the cleaning procedure can be started. No further contact to the agent is needed.



	Rinse	Cycles	Rinse Time	Toggle	Drying Time
	1 A	1 A	10 A	<input checked="" type="checkbox"/>	20 A
	2 A	1 A	10 A	<input checked="" type="checkbox"/>	240 A
+	A	A	A	<input type="checkbox"/>	A


Home Cancel Delete Row Previous OK

Volatile samples...

What if a sample is low-viscous and highly volatile?

Freely adjustable filling speed:

- For certain applications it is necessary to inject the sample slowly in the measuring system. By decreasing the fill rate filling errors can be prevented. The filling speed can be set easily.



Menu Meth... Meth... Dens... Xsample ... 2:42:59 PM Administrator

Xsample 340 Module Method Settings - Step 1 of 3 - General


Filling mode	Filling and Holding	
Fill volume	3.00	ml
Maximum fill rate	100	%
Fill volume for repetition	1.00	ml

Home Cancel Estimate Volume Next

Are there any special requirements to environmental conditions?

Developed for robustness:

- ▶ trouble-free operation are guaranteed due to the robust and streamlined design
 - ▶ Wetted parts made of high performance materials to ensure reliable operation
 - ▶ Valves proven in life science and chemical industry with highest chemical compatibility
- ... it can deal with almost everything.



Do I need to buy a new instrument when I want to automate my workflow?

Plug and play design:

- ▶ Connection to the master instrument via CAN cable is needed
- ▶ No additional housing for a module e.g. Lovis is needed any longer
- ▶ Small footprint as the unit can be placed beside or above the master instrument

... easy installation is garanted.



More about specifications.

FACTS & FIGURES

Specifications


Viscosity range	36.000 mPa.s at ambient and measuring temperature
Supported syringes	2 mL Luer plastic for sample viscosity up to 36,000 mPa.s 5 mL Luer plastic for sample viscosity up to 10,000 mPa.s 10 mL Luer plastic for sample viscosity up to 9,000 mPa.s 5 mL glass for sample viscosity up to 10,000 mPa.s
Filling mode	Pressure applied by syringe
Minimum sample volume per measurement	approx. 3 mL (with master instrument DMA M)
Sample recovery	not supported
Rinse consumption per cycle	approx. 7 mL (if rinse time: 10 s)
Compressed air supply (optional)	1.80–2.20 bar (26–32 psi) relative external compressed air typically: 2.0 bar (29 psi) relative

With or without cleaning?


XSAMPLE 340 – DIFFERENT MODELS

Two different models

- ▶ **Xsample 340**
 - no automated cleaning available
- ▶ **Xsample 340 with Cleaning**
 - Distribution plate
 - Cleaner unit




Xsample 340



- 1 Syringe adapter
- 2 Safety cover
- 3 Pusher drive


- 4 Distribution plate
- 5 Cleaner unit
- 6 Luer adapter



A look into the instrument's heart.

FILLING AND CLEANING PRINCIPLE

... with Cleaning




- ▶ Force controlled filling to ensure safe operation
- ▶ Connector for 2 bar external pressurized air
- ▶ Diaphragm pump for ambient air and cleaning agents
- ▶ Two cleaning agents

Xsample 340

SETTINGS

Methode settings:

- ▶ Filling mode
 - Filling and holding
 - Filling and releasing
- ▶ Fill volume
 - Estimate volume
- ▶ Maximum fill rate
- ▶ Fill volume for repetition




The screenshot shows the 'Xsample 340 Module Method Settings - Step 1 of 3 - General' window. The top menu bar includes 'Menu', 'Meth...', 'Meth...', 'Dens...', 'Xsample ...', and a status bar with the time '12:57:38 PM' and the user 'Administrator'. The main area contains four settings:

Setting	Value	Unit
Filling mode	Filling and Holding	
Fill volume	3.30	ml
Maximum fill rate	100	%
Fill volume for repetition	1.00	ml

At the bottom, there are five buttons: 'Home', 'Cancel', 'Estimate Volume', and 'Next'.

Filling mode – Filling and Holding:

- This mode keeps force on the plunger during measurement after sample has been filled.



Menu Meth... Meth... Dens... Xsample ... 12:57:38 PM Administrator


Xsample 340 Module Method Settings - Step 1 of 3 - General

Filling mode	Filling and Holding ▼		
Fill volume	3.30 A	ml	
Maximum fill rate	100 A	%	
Fill volume for repetition	1.00 A	ml	

Home Cancel Estimate Volume Next

Filling mode – Filling and Releasing:

- ▶ This filling mode releases the force on the plunger after sample has been filled.
This mode is recommended for high-viscous samples.



Menu Meth... Meth... Dens... Xsample ... 4:17:05 PM Administrator

Xsample 340 Module Method Settings - Step 1 of 3 - General


Filling mode	Filling and Releasing ▼		
Fill volume	3.30 A	ml	
Maximum fill rate	100 A	%	
Fill volume for repetition	1.00 A	ml	

Home Cancel Estimate Volume Next

Fill volume:

- ▶ Enter the volume of sample to be filled into your measuring system. The volume will be measured from the first contact of the pressure sensor to the syringe.

IMPORTANT: *Ensure that the entered filling volume corresponds with the volume of the connected measuring system. Otherwise the measuring system is not completely filled.*



Menu Meth... Meth... Dens... Xsample ... 12:57:38 PM Administrator

Xsample 340 Module Method Settings - Step 1 of 3 - General


Filling mode	Filling and Holding	
Fill volume	3.30 A	ml
Maximum fill rate	100 A	%
Fill volume for repetition	1.00 A	ml

Home Cancel Estimate Volume Next

Estimate Volume

- By tapping <Estimate Volume> the master instrument calculates automatically the minimum required sample volume for filling. This volume consists of the hose volume leading from Xsample 340 to the inlet of the master instrument and the volume of the measuring cell(s).

IMPORTANT: *The hose volume needs to be set to the correct value in the hardware configuration.*



Menu Meth... Meth... Dens... Xsample ... 12:57:38 PM Administrator


Xsample 340 Module Method Settings - Step 1 of 3 - General

Filling mode	Filling and Holding	
Fill volume	3.30 A	ml
Maximum fill rate	100 A	%
Fill volume for repetition	1.00 A	ml

Home Cancel Estimate Volume Next

Maximum fill rate:

- ▶ With this parameter you can set the rate filling volume per time. For certain applications it is necessary to inject the sample slowly in the measuring system. In this case decrease the fill rate.



Menu Meth... Meth... Dens... Xsample ... 12:57:38 PM Administrator

Xsample 340 Module Method Settings - Step 1 of 3 - General


Filling mode	Filling and Holding	
Fill volume	3.30 A	ml
Maximum fill rate	100 A	%
Fill volume for repetition	1.00 A	ml

Home Cancel Estimate Volume Next

Fill volume for repetition

- ▶ Xsample 340 performs the requested number of filling procedures by re-squeezing a specified sample volume through the measuring system without cleaning and drying in between.

IMPORTANT: *To avoid the syringe being emptied completely (which may lead to erroneous measurements), perform a test measurement .*




The screenshot shows the 'Xsample 340 Module Method Settings - Step 1 of 3 - General' window. The 'Fill volume for repetition' field is highlighted with a red box, showing a value of 1.00 A ml. Other fields include 'Filling mode' (Filling and Holding), 'Fill volume' (3.30 A ml), and 'Maximum fill rate' (100 A %).

Setting	Value	Unit
Filling mode	Filling and Holding	
Fill volume	3.30 A	ml
Maximum fill rate	100 A	%
Fill volume for repetition	1.00 A	ml

Methode settings:

- ▶ Draining mode
 - With Air
 - With External Compressed Air
- ▶ Draining Time


It is essential that the sample is drained out completely, before the cleaning procedure with rinsing agent is started.



The screenshot shows the 'Xsample 340 Module Method Settings - Step 2 of 3 - General' window. The interface includes a menu bar at the top with options: Menu, Meth..., Meth..., Dens..., Xsample ..., and a status bar showing the time 2:17:22 PM and the user Administrator. The main content area contains three settings: 'Auto air check before filling' with an unchecked checkbox, 'Draining mode' set to 'With Air' in a dropdown menu, and 'Draining Time' set to '30' in a text field with a unit selector 'A' and a small 's' to the right. At the bottom, there are five buttons: Home, Cancel, a disabled button, Previous, and Next.

Cleaning settings:

- ▶ Toggle:
 - Check box is activated:
During the rinsing procedure small air-liquid sections (air gaps) will be generated in the hose to ensure better cleaning performance.
 - Check box is deactivated:
Xsample 340 is rinsing the connected measuring system without air-liquid sections.



	Rinse	Cycles	Rinse Time	Toggle	Drying Time
1	A	1 A	10 A	<input checked="" type="checkbox"/>	20 A
2	A	1 A	10 A	<input checked="" type="checkbox"/>	240 A
+	A	A	A	<input type="checkbox"/>	A

Home Cancel Delete Row Previous OK

Ensure that before cleaning procedure is started, the measuring system is free of sample!

THANKS FOR YOUR ATTENTION