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

USP Dietary Supplements Compendium 2015

Two Volume set. One indispensable quality.
Order the updated and expanded 2015 Dietary Supplements

FEATURING

- 75 new dietary supplement monographs--nearly 500 in all--from USP 38-NF 33 through the First Supplement.
- 27 new General Chapters
- More than 175 excipient monographs
- Over 200 Food Chemicals Codex (FCC) monographs
- More than 40 new and revised DSC admission evaluations
- Includes over 150 added pages of color plates and illustrations

DOWNLOADS

- 2015 DSC Volume 1: Table of Contents and Reference Standards Index
- 2015 DSC Volume 2: Tables of Contents

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7044 (203) High-Performance Thin-Layer Chromatography / Chemical Tests First Supplement to USP 38–NF 33

Add the following:

•(203) HIGH-PERFORMANCE THIN-LAYER CHROMATOGRAPHY PROCEDURE FOR IDENTIFICATION OF ARTICLES OF BOTANICAL ORIGIN

INTRODUCTION

This chapter describes a procedure for use in a USP identification test that relies on the technique of high-performance thin-layer chromatography (HPTLC). It is applicable to the identification of articles of botanical origin in USP's compendia that serve as a drug, substance or drug product, or as an ingredient or a dietary supplement. Careful control of the variables for the HPTLC technique is briefly described with references to more detailed information provided in equipment manuals. Reproducibility of the results allows comparison of closely related botanical materials that are not official ingredients. The analytical technique uses high-performance plates, appropriate equipment to control variables, and a system suitability test for purposes of performance qualification. [NOTE—See Identification of Articles of Botanical Origin Using High-Performance Thin-Layer Chromatography Procedure (1064) for additional information.]

REQUIRED EQUIPMENT

The equipment used for HPTLC technique consists of the following:

- Plates: Unless otherwise specified in the individual monograph, use plates coated with a uniform 200-µm layer of porous (10-Å pore size) silica gel having regular particles of 2–10 µm and an average particle size of 5 µm, a polymer binder, and a fluorescence indicator (F_{254} 254 nm). [NOTE—Chromatographic methods using high-performance thin-layer chromatographic glass plates are preferred over aluminum-backed sheets because of greater mechanical stability.]
- A device suitable for the application of specified volumes of samples as bands with specified length at the specified positions.
- A suitable chromatographic chamber (for example, a twin trough chamber) allowing for control of saturation and developing distance.
- A device suitable for controlling the activity of the stationary phase via relative humidity.
- A device suitable for reproducible drying of the developed plate.
- A device suitable for treatment of the plate with derivatization reagent, if required.
- A device suitable for heating as part of the derivatization procedure, if required.
- A system suitable for documentation of chromatograms under UV 254 nm, UV 366 nm, and white light.

Each of these devices as well as the system as a whole should pass installation qualification (IQ), operational qualification (OQ), and performance qualification (PQ) in order to ensure that the instruments are working according to their specifications to control the variables within their intended ranges. [NOTE—See Analytical Instrument Qualification (1058) for additional information.] Installation qualification and operational qualification are typically performed by the instrument manufacturer. The system suitability test is performed by the analyst as proof of performance qualification.

PROCEDURE

Preparation of Test Solution

Unless otherwise stated in an individual monograph, 100 mg of a powdered botanical ingredient, 10 mg of a dry extract or fraction, or the amount of a dosage form containing the equivalent of the aforementioned quantity of the botanical ingredient is extracted for 15 min with 1 mL of methanol. After centrifugation the filtrate or supernatant is used as the sample solution. Unless otherwise stated in an individual monograph, 10 µL of an essential oil is dissolved in 1 mL of toluene and used as the sample solution.

Preparation of the Standard Solutions

Unless otherwise stated in an individual monograph, USP Reference Standards of individual marker compounds are dissolved at a concentration of 1 mg/mL in methanol. The USP Reference Standard extracts are shaken and sonicated in methanol at a concentration of 10 mg/mL; for essential oils, the USP Reference Materials are dissolved in toluene at a concentration of 10 µL/mL.

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b. Thin Layer Chromatography


Fig. 5 Typical HPTLC chromatograms

Track assignment: 1) USP Piperine RS (0.9 mg/mL), Borneol (1.0 mg/mL), with increasing R_f ; 2) USP Black Pepper Powdered Extract RS (25.0 mg/mL); 3–5) Black Pepper (commercial samples); 6) Black Pepper (authenticated sample)

Sample solutions:	according to the monograph
Standard solutions:	in methanol
Plate:	HPTLC, Si 60 F ₂₅₄
Developing solvent:	hexanes and ethyl acetate (5:3)
Developing distance:	6 cm
Saturation time:	20 min, with filter paper
Relative humidity:	33%, saturated MgCl ₂
Temperature:	ambient, not to exceed 30°
Application volume:	7 µL of the Sample solution and 3 µL of the Standard solutions, as 8-mm bands
Detection:	a) underivatized, under UV light at 254 nm b) derivatized, under white light
Derivatization reagent:	<i>p</i> -anisaldehyde, glacial acetic acid, methanol, and sulfuric acid (0.5:10: 85:5)

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

Results:

Fig. 1 Calibration function of metronidazole (X) and determination in drug sample (+)
(P200-04 090814-03)

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



First application position X: 30.0 mm

Application position Y: 8.0 mm

Distance between tracks: 14.0 mm ☒ Automatic ☐ Manual



Track	Application position	Application volume	Rack col	Rack row
1	30.0 mm	2.0 μ L	A	1
2	44.0 mm	3.0 μ L	A	2
3	58.0 mm	3.0 μ L	A	2
4	72.0 mm	4.0 μ L	A	1
5	86.0 mm	5.0 μ L	A	3
6	100.0 mm	5.0 μ L	A	3
7	114.0 mm	6.0 μ L	A	1
8	128.0 mm	7.0 μ L	A	4
9	142.0 mm	7.0 μ L	A	4
10	156.0 mm	8.0 μ L	A	1
11	170.0 mm	10.0 μ L	A	1

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
Fig. 1 60 samples (6 batches, B1-B6) and 12 standards* are applied onto one HPTLC plate; Coenzyme Q10 is seen at $R_f \approx 0.20$. In the Q10-batches on the center tracks vitamin B2 and vitamin E are also present ($R_f \approx 0$ and 0.4); (NOTE: For better visualization in this figure the application volume was increased to $6 \mu\text{L}$).

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



Densitogram of eight phospholipids

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

Figure 1. Offish, *Ruvettus pretiosus*

Figure 2. Escolar, *Lepidionotus flavobrunneum*

Figure 3. Rockfish, *Centrolophus niger*


Figure 4. Rockfish, *Tudora sp.*

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
▲ Selected tracks are detected under UV 365/>400 nm: Track 1 (standard mixture of the 7 biogenic amines β -phenylethylamine, tyramine, histamine, spermine, spermidine, cadaverine and putrescine); Track 2 (not contaminated fish sample spiked with standard mixture); Track 3 (as Track 2, but developed with solvent proportions 10:6:2).

► S. Kretzschmar and K. Speer in CBS 95

▲ Densitogram of a sample "Tuna in oil" (1:10 diluted) at UV 365/>400 nm

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
The screenshot shows the VisionCATS software interface. The top navigation bar includes tabs for HPTLC-Steps, Chromatography, Data, Spectrum, SST, and Evaluation (6.1). The Evaluation Steps section on the left shows a workflow from Definition to Results, with the Results step highlighted by a red circle. The main Results section displays data for Sinensetin @ 366 nm, with 8 samples assigned. The data is organized into four groups, each representing a sample, with sub-rows for Volume and individual tracks.

Sinensetin @ 366 nm (8 samples assigned)				
Sample 'R3518_150106_02'				
		102.3 ng/ml	CV=0.365 %	(2 applications) 102.3 µg in 1.00 mg
Volume: 5.0 µl		102.3 ng/ml	CV=0.365 %	(2 replicas)
Track 7	Rf 0.331	102.6 ng/ml	513.0 pg	
Track 15	Rf 0.330	102.1 ng/ml	510.3 pg	
Sample 'S10795_150106_01'				
		97.66 ng/ml	CV=0.824 %	(2 applications) 48.83 µg in 500.00 mg
Volume: 5.0 µl		97.66 ng/ml	CV=0.824 %	(2 replicas)
Track 1	Rf 0.313	98.23 ng/ml	491.2 pg	
Track 9	Rf 0.334	97.10 ng/ml	485.5 pg	
Sample 'S10795_150106_02'				
		71.61 ng/ml	CV=0.339 %	(2 applications) 35.80 µg in 500.00 mg
Volume: 5.0 µl		71.61 ng/ml	CV=0.339 %	(2 replicas)
Track 3	Rf 0.317	71.78 ng/ml	358.9 pg	
Track 11	Rf 0.334	71.44 ng/ml	357.2 pg	
Sample 'S10795_150106_03'				
		62.51 ng/ml	CV=0.155 %	(2 applications) 31.26 µg in 500.00 mg
Volume: 5.0 µl		62.51 ng/ml	CV=0.155 %	(2 replicas)
Track 5	Rf 0.325	62.45 ng/ml	312.2 pg	
Track 13	Rf 0.333	62.58 ng/ml	312.9 pg	

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The screenshot shows the CAMAG TLC Visualizer software interface. The top window, titled '5.1 Evaluation', is highlighted with a yellow circle. It displays a list of substances and their corresponding peak data. The interface includes a menu bar with options like 'TLC-Steps', 'Chromatography', 'Data', 'SST', and '5.1 Evaluation'. The main content area shows a table of results for two substances: 6-Gingerol and 6-Shogaol. Each substance has two sample assignments with their respective concentrations and standard deviations.

Substance	Sample	Concentration (µg/ml)	σ (%)	Applications
6-Gingerol	Sample 'S13247'	662.264	0.000	1
	Sample 'S13248'	741.627	0.000	1
6-Shogaol	Sample 'S13247'	204.401	0.000	1
	Sample 'S13248'	329.587	0.000	1

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

The screenshot displays the visionCATS software interface. The main window shows a 'Sequence' table with columns 'Tr.', 'Description', and 'Vol.'. A 'TLC Steps' panel is visible, showing a sequence of steps with icons for sample application, development, and detection. A 'Sequence Table Editor' dialog box is open, showing a list of samples with their 'Tr.' and 'Vial ID'. A 'Scanner Settings' dialog box is also open, showing configuration options for the scanner.

Sequence Table Editor:

Tr.	Vial ID
1	R13360-01-9%
2	R13360-01-16%
3	R13360-01-33%
4	R13360-01-50%
5	R13360-01-100%
6	R13359-01-9%
7	R13359-01-16%
8	R13359-01-33%
9	R13359-01-50%
10	R13359-01-100%
11	S13247
12	S13248
13	S13249
14	S13250
15	S8122
16	

Scanner Settings:

- Scanner type: Single λ
- Optimization for: Resolution
- Measurement mode: Absorption
- Speed/resolution/slit: User
- Wavelength selection:
 - Selection: User
 - Wavelength(s) applied: 232 nm
 - Lamp: Deuterium & Tungsten

The 'Wavelength selection' section includes a graph showing the spectral intensity of the Deuterium and Tungsten lamps. The x-axis represents wavelength in nm (200 to 900), and the y-axis represents intensity. The Deuterium lamp shows a strong peak around 232 nm, while the Tungsten lamp shows a broad peak around 600 nm.

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Ginger rhizome, Jiang (*Zingiber officinale*)

Related documents:

SOP 001 General Methodology for HPTLC
SOP 002 Identification by HPTLC

1. Scope
This method identifies dried Ginger rhizome (*Zingiber officinale* Roscoe) by HPTLC fingerprint and discriminates dried Lesser galangal rhizome (*Alpinia officinarum* Hasske), dried rhizome of *Kaempferia galanga* (*Kaempferia galanga* L.), and dried *Kaempferia galanga* (*Kaempferia galanga* L.) rhizome.

2. Source of method
Monograph Ginger Rhizome, USP Documentation

3. Procedure
Sample preparation:

Reference substances:

Tr.	Description	Vol.
1	6-Gingerol	2.0
2	6-Shogaol	2.0
3		
4		
5		
6		
7		
8		
9		
10		
11		
12		
13		
14		
15		

Stationary phase:

Application:

Mobile phase:

Development:

Derivatization reagent:

Documentation:

Created by: EC/18 Nov 2009

4. Results
Note: The images presented in this section are examples and are not intended to be used as basis for setting specifications for quality control purposes.

Fig. 1) *Alkaldehyde* reagent, white RT

Identification
Compare result with reference images. The fingerprint of the test solution is similar to that of the corresponding botanical reference sample. Additional weak zones may be present.

Test for adulteration
Under UV 366 nm no green or blue fluorescent zones are seen between Rf ~ 0.2 and 0.4 (Lesser galangal rhizome). No green fluorescent zone is seen just below the position of reference substance 6-gingerol (Rf ~ 0.2). No green fluorescent zone is seen

Method: Ginger rhizome

Method: Ginger rhizome

Zingiber officinale

40 distinct methods on CAMAG Server 44 methods on CAMAG Server 6 distinct methods downloaded 8 methods downloaded

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

Request for Laboratory Services

Information about customer	
Name of company	
Person to contact	
Responsibility within the company	
Department	
Address	
Phone	
Fax	
E-Mail	
Field of activity/business	
Information about analytical task (HPTLC)	
Type of Service requested:	
☒	feasibility study
☐	method development
☐	validation
Detailed description of the analytical goal:	
Current method of analysis (please attach method):	
☒	TLC/HPTLC
☐	HPLC
☐	GC
☐	Other
☐	None
Type of sample (matrix):	
☒	pharmaceutical formulation (tablet, syrup, etc.)
☐	cosmetic
☐	plant extract
☐	oil/glycerin
☐	other (specify):

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[TLC / HPTLC](#)
Thin-Layer Chromatography

[DBS](#)
Dried Blood Spot Sampling

[ALOX](#)
Aluminium Oxide

[CAMAG](#)
About us

INSTRUMENTS, TOOLS AND CONCEPTS FOR HPTLC

CAMAG is the world leader in instrumental Thin-Layer Chromatography (TLC/HPTLC).

We provide high-end instruments, tools and concepts for TLC/HPTLC. Additionally, we support our clients with training, technical and contract laboratory services.

Wherever complex substance mixtures are analyzed (e.g. herbal, environmental or forensic samples), TLC/HPTLC often is an excellent alternative to GC and HPLC. Moreover, it has proven value as a solid quantitative analysis technique in quality assurance of food, pharmaceuticals and other industrial products.

Since 1961, CAMAG is dedicated to the development and manufacturing of instruments, software and consumables for all steps of the TLC/HPTLC procedure. Our products are truly Swiss made and have an excellent reputation throughout the world.

NEWS

19.11.2015 - New case study now online!
A case study from the field of forensic drug control is now online: >> In-process control during chemical synthesis...

UPCOMING EVENTS

PITTCON 2016
06. - 10. March 2016
Atlanta, GA (USA)

ARABLAB 2016
20. - 23. March 2016
Dubai, United Arab Emirates

[+ more](#)

METHODS / APPLICATION NOTES

HPTLC Identification of Fatty ...
Determination of aflatoxins B1...

Ginger
Feverfew

[+ more](#)

LATEST PUBLICATIONS

Journal of Ethnopharmacology: ...
Food Additives & Contamina...
Fast UHPLC-MS/MS method to sim...
THIN-LAYER CHROMATOGRAPHY | Ov...

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