

## FILMDIST® SP 200



[www.pilodist.de](http://www.pilodist.de)

Tel: ++ 49 (0) 2225 955910

Fax: ++ 49 (0) 2225 9559111

e-mail : [info@pilodist.de](mailto:info@pilodist.de)

PILODIST GmbH  
Eichelnkampstraße 2  
D-53340 Meckenheim  
Sitz: Bonn

Eingetragen AG Bonn HRB 12941  
Managing Director Rolf Jeub  
USt.-ID.-Nr. DE 813 408 398

Sparkasse KölnBonn  
Kto. 33302365  
BLZ 37050198  
IBAN DE 05 3705 0198 0033 3023 65  
Swift-BIC: COLSDE33

## FILMDIST® SP 200

Universal lab-scale apparatus for shortpath evaporation.

Rotary thinfilm evaporator with thermostat mantle and internal condenser for continuous shortpath distillation or gentle concentration of thermally sensible products. Suited for laboratory or pilot operation. The pilot operation data and experimental values can be transferred to bigger systems.

The standard system uses glass wiper elements. Other wiper elements are available.

### *Technical Data*

Operating temperature:	max. 250°C
Operating pressure:	0.001 - 1000 mbar
Feed range:	50 - 1000 ml / h
Evaporator surface:	200 cm <sup>2</sup> , short path
Evaporator material:	borosilicate glass
Feed vessel:	250 ml
Product receivers:	250 ml
Mains supply:	208 - 250 V, 50 Hz (60 Hz upon request)

[www.pilodist.de](http://www.pilodist.de)

Tel: ++ 49 (0) 2225 955910

Fax: ++ 49 (0) 2225 9559111

e-mail : [info@pilodist.de](mailto:info@pilodist.de)

PILODIST GmbH  
Eichelnkampstraße 2  
D-53340 Meckenheim  
Sitz: Bonn

Eingetragen AG Bonn HRB 12941  
Managing Director Rolf Jeub  
USt.-ID.-Nr. DE 813 408 398

Sparkasse KölnBonn  
Kto. 33302365  
BLZ 37050198  
IBAN DE 05 3705 0198 0033 3023 65  
Swift-BIC: COLSDE33